
A 216 cm superabsorbent floormat - easy to fold and form

DryMax Floor mat TRIPLE



Place it in front of or beside the feet of the surgical team. During absorption, the mat swells and its contents turn to gel.

The mat is dual-sided and can be placed on floors that are already wet. If required, additional mats can be added or the product can be used in conjunction with DryMax Floor mat 2.4 to provide extra absorption.

Proper fluid management promotes Health and Safety as it reduces the risk of slips and trips and provides a cleaner and safer working environment.

As DryMax Floor mat TRIPLE locks in the contaminated fluids without dripping, it helps reduce the risk of spreading contaminants. In this way, the mat contributes to hygiene control and staff safety.

Cleaning times and costs are also significantly reduced as fluids are absorbed and retained by the mat.

ltem	REF	Pcs/box	Product size	Absorbency		
				Tap water	NaCl 0,9%	Packaging
DryMax Floor mat TRIPLE	F20243	40	216 cm x 37 cm	15 l	61	40 pcs/box













Absorbs water-based liquids such as saline, urine, amniotic fluids and blood. In the case of liquids containing bodily fluids and proteins, the absorption capacity will decrease compared to the saline data provided.





- **Operating Rooms**
- **Examination Rooms**
- Gynaecology
- Urology
- **Emergency rooms**
- **Ambulances**
- Body bags
- Staff Safety
- **Hygiene Control**
- Faster turnover time
- Easy to form and fold

